
CCL VARIO SYSTEM

Corneal Crosslinking
3 mW – 9 mW – 18 mW

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Integriertes Display Integrated monitor

Transportkoffer Sturdy transport case

CCL vario

CCL – Corneal Crosslinking

Corneal crosslinking (CCL) is a treatment with the aim of strengthening the corneal stromal tissue through the formation of new chemical bonds between stromal fibers.

CCL is the only effective treatment to stop progressive keratoconus as well as related ectatic disorders (such as pellucid marginal degeneration and iatrogenic ectasia) and has a regularization effect on corneal topography.

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Clinical experience

In recent years corneal crosslinking has become the standard procedure for treating patients with progressive keratoconus and other ectatic corneal diseases because of its effectiveness and lack of serious side effects.

In a number of clinical studies it was demonstrated that in more than 85 % of eyes treated the BCVA increased significantly. Six months after the procedure cylinder was reduced in the majority of all patients. The steepest K-value was reduced by an average of up to two dioptries.

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
Background

Corneal crosslinking is a process of photopolymerization. During this process singlet oxygen is being created with the use of riboflavin as a photomediator activated by UV-light. Free radicals lead to physical intra- and interhelical crosslinks of stromal collagen fibers.

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
CCL vario – device details

CCL vario comes with one diode and special optics which is homogenizing the beam. Thus hot spots are being avoided and the endothelium is sufficiently protected. To offer more flexibility, it enables the surgeon to choose between three energy levels. This allows the surgeon to choose the ideal energy/time combination for the intended treatment. At the end of the procedure the system switches off automatically. Since the total dose of radiation is time dependent, each setting delivers a total energy of 5,400 mJ.

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- three energy levels allow high flexibility
- only five minutes treatment time in the highest energy level
- continuously adjustable aperture for protection of limbal stem cells
- efficient and safe to use
- user friendly
- sturdy transport case allowing the surgeon to be mobile

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Wavelength range:	365 nm
Illumination intensity:	3, 9 or 18 mW/cm ²
Working distance:	50 mm ± 5 mm
Light emission:	Continuous wave (CW)
Spot sizes (continuously adjustable):	7–11 mm
Timer:	30, 10 or 5 min
Electric power:	100–240 V
Dimensions Hard case (W x L x H):	37 cm x 46 cm x 14 cm
Weight (total):	6.2 kg


CE 1275

Technical data

To guarantee the high level of safety the beam of the CCL vario has a waistline at a distance of 50 mm from the optics and a depth of focus of approximately +/- 5 mm.

To protect the limbal stem cells and to focus the beam on the clear cornea the CCL vario has a continuously adjustable aperture from 7 mm to 11 mm.

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Reliability, the highest quality, service-oriented actions and enthusiasm for our work – for more than 15 years, this is what AIVIMED has stood for.

AIVIMED GmbH
medical optics

Borsigstr. 34
65205 Wiesbaden
Germany

T +49 (0) 611 450 32 0
F +49 (0) 611 450 32 20

info@aivimed.com
www.aivimed.de
www.aivimed.com